

Pacific Winemaking LLC

Tannin Bench Trials

1. Make up a batch of 90ml Distilled H₂O to 10ml Ethanol for each tannin sample to be trialed. This is A.
2. Mix in 1 gram of the Tannin into 100ml of A. This creates a 1% solution.
3. Set up a minimum of 6 tasting glasses per sample for trial. Mark one glass as a control and ensure no additions of tannin are made to this glass.
4. Add 100ml of the wine to be trialed to each of the 6 tasting glasses.
5. Based on 1% Tannin solution added to 100ml of wine, 0.1ml is equivalent to 1 g/hL addition of the Tannin.
6. Utilize the Tannin Product Data Sheet to determine the suggested range of additions for the wine to be treated.
7. Add 5 different levels of Tannin to provide a comparison of effect.

Equipment Needed:

6 x tasting glasses
1 x 600ml Beaker
1 x 1ml Graduated Pipette
Distilled Water
Ethanol
Tannin sample(s)

Please contact us for larger trials and recommendations of suitable products.

PO BOX 731, CARLTON, OREGON 97111

Ph: +1(503) 419-7942 **Fax:** +1(510) 764-9713 info@pacificwinemaking.com

www.pacificwinemaking.com